

Fiber coupled FTIR-spectroscopy for biomedical diagnostics

V. Artyushenko*, J. Mannhardt*, N. Pleshko**, C. McGoverin**, Q. Onigbanjo**, T. Sakharova***
 *art photonics GmbH, Rudower Chaussee 46, 12489 Berlin, Germany, sa@artphotonics.com
 **Dept. of Bioengineering, Temple University, Philadelphia, PA, USA
 *** General Physics Institute RAN, Moscow, Russia



Development of Mid IR fiber optics within the last decade enables biomedical applications of FTIR-spectroscopy not only for *in-vitro*, but for *in-vivo* diagnostics. While in the past fiber coupled FTIR-spectrometers were produced with LN-cooled MCT detectors, now very good S/N Ratio can be reached with room temperature DTGS-detectors. High quality Polycrystalline PIR-fibers from Silver Halide crystals and Chalcogenide As-S-glass CIR-fibers provide good throughput in optimal coupling with FTIR - making biomedical diagnostics possible using ATR absorption and reflectance spectroscopy in Mid IR - the most informative finger-print region for organic molecular vibrations.



FTIR spectrometers equipped by DTGS-detectors & coupled with ATR-fiber probes by mirror-fiber couplers installed in sample chambers.
 Fig. 1. Fiber coupled FTIR from Thermo (iS5) and ABB

While all fiber optic probes reduce S/N Ratio when coupled to FTIR-spectrometers - due to an evident troubles of efficient coupling, - their spectra quality is almost the same as for a common ATR-accessories used for FTIR with DTGS detectors:




Fig. 2. Example of Acetone spectra measured with iS5 FTIR from Thermo equipped with standard ATR-accessories and 1,5m long PIR-fiber ATR-Probe

At the 1st stage fiber coupled FTIR-spectrometers should be used for intensive studies of tissue or bioliquid in a broad spectral range to detect specific changes related to the disease, but these changes are already known for the set of few wavelengths (signal & reference) - then small and low cost fiber sensors can be developed for specific diagnostics.




Fig.3. Fabry-Perot Spectral Engine from VTT vs iS5

ATR-Probes


ATR-Absorption Spectroscopy with detachable PIR-fiber loops enable molecular analysis of liquids and tissue *in-vivo*, in real-time, including endoscopic applications, in two parts of Mid IR-range (due to use of CIR- or PIR-fibers): 6500-1700cm⁻¹ and 3600-600cm⁻¹.

Fig. 4. Transmission spectra of different ATR-probes of 1,5m length and design of detachable single & multi-loop ATR-tips.




ATR-probes with detachable tips based on PIR-Loops (or other ATR-tips - like Silicon or Diamond-cone) enables to sterilize them and to use for diagnostics *in-vivo* (which is also needed for disposable PIR-loops).

Fig. 5 Examples of ATR-absorption spectra for a hand skin.




Sterilization is needed for all biomedical applications *in-vivo* and it was tested for detachable PIR-Loop ATR-tips. Spectral overlay of bovine cartilage after sterilization does not show any significant changes in the Amide I & II peak heights

Fig. 6. Photos of PIR-loops glued in detachable PEEK-caps & Absorption spectra of bovine cartilage after PIR-Loops sterilization



FTIR-spectroscopy with ATR-Fiber Probes in biomedical applications

Pioneering applications of FTIR-fiber spectroscopy were started long time ago to detect malignant tissue during cancer operation (#1) and now much more trials are done *in-vivo* (#2) & *in-vitro* (#3) - see data presented from 3 articles at Fig. 7 to show great potential of FTIR-spectroscopy for molecular tissue diagnostics.



#3 1 October 2010 / Vol. 1, No. 3 / BIOMEDICAL OPTICS EXPRESS 1014





Fig. 8. Representative FTIR spectra for EUSC study and normal (black) colonic mucosa over fingerprint region (900 to 1800 cm⁻¹) from the right (R), transverse (T), and left (L) anatomic segments demonstrate very good SNR and several numerous absorbance peaks in 6 separate sub-ranges.



#3. Spectra of normal blue, necrotic (green) and carcinoma (red) brain tissue. N.Bergner et al., Analyst, 2013, 138, 3983